



## > Catalase Fluorometric Detection Kit

Catalog # 907-027  
Sufficient Reagents for 500 Tests



Store at 4°C



Check our website for additional protocols, technical notes and FAQs



For proper performance, use the insert provided with each individual kit received

### Table of Contents

<u>2</u>	<u>Introduction</u>
<u>3</u>	<u>Materials Supplied</u>
<u>3</u>	<u>Storage</u>
<u>3</u>	<u>Materials Needed but Not Supplied</u>
<u>4</u>	<u>Reagent Preparation</u>
<u>5</u>	<u>Sample Handling</u>
<u>5</u>	<u>Assay Procedure</u>
<u>6</u>	<u>Calculation of Results</u>
<u>6</u>	<u>Typical Results</u>
<u>7</u>	<u>References</u>
<u>8</u>	<u>Limited Warranty</u>

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

## Introduction

Assay Designs' Catalase Fluorometric Detection Kit is a sensitive fluorescent assay to detect catalase activity by measuring the amount of substrate (hydrogen peroxide) remaining after sample addition. This kit provides a simple homogenous assay that is adaptable to kinetic and high throughput applications.


Catalase is an antioxidant enzyme that catalyzes the decomposition of hydrogen peroxide ( $H_2O_2$ ) to water and oxygen. Catalase is ubiquitously expressed in mammalian and non-mammalian aerobic cells containing the cytochrome system and is a protein tetramer of four ferrihemoprotein groups per molecule. The enzyme has been isolated from various sources, including bacteria and plant cells<sup>1-3</sup>. Catalase activity varies greatly from tissue to tissue. Highest activity is seen in liver and kidney, while lowest activity is seen in connective tissue<sup>3</sup>. The presence of catalase in the peroxisomes of mammalian cells helps to alleviate oxidative damage by catalyzing the metabolism of the peroxide substrate<sup>4</sup>.

The production of  $H_2O_2$ , a highly reactive byproduct of mammalian metabolism, can be deleterious to cells and is an end product result of various oxidase and superoxide dismutase reactions. Accumulation of  $H_2O_2$  can result in cellular damage through oxidation of proteins, DNA, and lipids thus resulting in cell death and mutagenesis<sup>8-11</sup>. The role of  $H_2O_2$  in oxidative stress related pathologies (e.g., inflammation, cancer, diabetes, cardiovascular disease, anemia, Parkinson's disease, Alzheimer's disease) has been widely studied<sup>8,12</sup>.

## Principle

The Catalase Fluorometric Detection Kit is a sensitive assay that utilizes a non-fluorescent detection reagent to measure  $H_2O_2$  substrate left over from the catalase reaction<sup>5-6</sup>.

Reaction:



## **Materials Supplied**

- 1. Catalase Reaction Buffer Concentrate**  
20 mL, Product No. 80-1697  
A 5X Buffer Concentrate
- 2. Detection Reagent**  
1 vial, Product No. 80-1698  
One vial is sufficient for 500 tests
- 3. Horseradish Peroxidase Concentrate**  
1 vial, Product No. 80-1699  
18.9 units of enzyme
- 4. Hydrogen Peroxide**  
200  $\mu$ L, Product No. 80-1700  
A 3% stabilized solution of hydrogen peroxide
- 5. Catalase Enzyme**  
1 vial, Product No. 80-1701  
A crystalline suspension of active catalase enzyme

## **Storage**

All components of this kit are stable at 4°C.

## **Materials Needed but Not Supplied**

1. Dimethyl sulfoxide (DMSO)
2. Black 96 well plates
3. Fluorescence plate reader
4. Deionized water
5. Superoxide dismutase (optional, see Reagent Preparation #3)

## Reagent Preparation

### 1. 1X Reaction Buffer

Prepare a 1X solution of the Reaction Buffer by diluting 4 mL of Reaction Buffer Concentrate with 16 mL of deionized water. This can be stored at 4°C until the kit expiration date, or for 3 months, whichever is earlier.

### 2. 100X HRP

Prepare a 100X solution of Horseradish Peroxidase (HRP) by diluting it 1:94.5 with 1X Reaction Buffer. For example, add 5  $\mu\text{L}$  of the HRP stock reagent provided to 467.5  $\mu\text{L}$  1X Reaction Buffer. Make enough HRP for one day worth of experiments. This 100X stock may be stored for several weeks at or below -20°C. Aliquot in single use vials.

### 3. Detection reagent

Dissolve the contents of one vial in 500  $\mu\text{L}$  DMSO. Allow the contents to sit at room temperature for 15 minutes. Gently pipet up and down several times to dissolve any clumps. Once dissolved the detection reagent should be used promptly and any remaining reagent can be aliquotted and stored at -70°C. Avoid repeated freeze thaw cycles.

**At NADH levels above 10 $\mu\text{M}$  and Glutathione levels above 300 $\mu\text{M}$ , detection reagent oxidation results from side chain reaction between NADH and HRP. This could result in aberrant readings. To minimize this interference it is recommended to add superoxide dismutase (SOD) at 40U/mL to the reaction<sup>7</sup>.**

### 4. Hydrogen Peroxide (40 $\mu\text{M}$ )

The 3% hydrogen peroxide solution provided is ~ 0.881 M. Prepare a 40  $\mu\text{M}$  solution of the hydrogen peroxide in 1X Reaction Buffer. This will serve as the substrate for catalase.

### 5. Catalase Standard Curve

See vial for specific catalase activity and concentration. The Catalase is a crystalline suspension in water. This crystalline suspension must be dissolved before use. Vortex the vial to evenly disperse the crystalline suspension. Pipet an appropriate amount of catalase into 1X Reaction Buffer. Warming gently at 30°C and slight agitation will help dissolve the catalase crystals.

One unit will decompose 1.0  $\mu\text{mole}$  of  $\text{H}_2\text{O}_2$  per min at pH 7.0 at 25°C, while the  $\text{H}_2\text{O}_2$  concentration falls from 10.3 to 9.2 mM, measured by the rate of decrease of  $A_{240}$ .

Prepare a standard curve by diluting the appropriate amount of Catalase in 1X Reaction Buffer. A suggested standard curve range is 0 to 4 units/mL.



The fluorescent product of the detection reagent is not stable in the presence of thiols (DTT or 2-mercaptoethanol). Keep these reactants below 5  $\mu\text{M}$ .

## 6. Reaction Cocktail

Prepare 10 mL **Reaction Cocktail** for 100 assays:

- 100 µL Detection Reagent .
- 100 µL 100X HRP.
- 9.8 mL of 1X Reaction buffer

The volume of the **Reaction Cocktail** can be scaled down or up (as needed) provided that the ratios of the ingredients are kept constant.

The fluorescent product of the detection reagent is not stable in the presence of thiols (e.g. DTT or 2-mercaptoethanol). Keep these reactants below 5 µM. If using your own buffer, keep the reaction between pH 7-8 (pH 7.4 is optimal).

## Sample Handling

Dilute samples in 1X Reaction Buffer. If the general amount of catalase activity is not known, make several dilutions (1, 10, 20, 50 fold dilution) of your sample. One of these dilutions should fall in the standard curve range.

## Assay Procedure

1. Pipet 50 µL of standard or sample to the bottom of a 96 well black plate.
2. Pipet 50 µL of the 40 µM H<sub>2</sub>O<sub>2</sub> solution to each well.
3. Incubate the plate for 30 to 60 minutes at room temperature.
4. Add 100 µL of the Reaction Cocktail to each well.

**Note:** Each investigator should optimize incubation times for their particular application.

5. Incubate the plate for 10 to 15 minutes.
6. Read samples using excitation 530-570 nm (570 nm is optimal) and measure fluorescence at 590-600 nm.



All standards, controls, and samples should be run in duplicate.



Pre-rinse each pipet tip with reagent. Use fresh pipet tips for each sample, standard, and reagent.



Pipet the reagents to the sides of the wells to avoid possible contamination.



Make sure to multiply sample concentrations by the dilution factor used during sample preparation.

## Calculation of Results



Several options are available for the calculation of Catalase activity in samples. We recommend that the data be handled by an immunoassay software package utilizing a 4 parameter logistic curve fitting program. If data reduction software is not readily available, the Catalase activity can be calculated as follows:

1. Calculate the average net RFU for each standard and sample by subtracting the average 0 U/mL RFU from the average RFU for each standard and sample.
2. Using linear graph paper, plot the Average Net RFU for each standard versus Catalase Activity (U/mL) in each standard. Approximate a line through the points. The concentration of the unknowns can be determined by interpolation.

Samples with concentrations outside of the standard curve range will need to be reanalyzed using a higher dilution.

## Typical Results

The results shown below are for illustration only and should not be used to calculate results from another assay. The reaction contained 20  $\mu\text{M}$   $\text{H}_2\text{O}_2$  (final) per well and the indicated amounts of catalase in 1X reaction Buffer. The reaction was incubated for 30 minutes at room temperature. The graph reports the change in fluorescence, observed fluorescence from negative control (no catalase) minus catalase sample fluorescence.



## References

1. Deisseroth, A., and Dounce, A.L. (1970) *Physiol Rev.* 50, 319-375.
2. Mueller, S., Riedel, H. D., and Stremmel, W. (1997) *Anal Biochem.* 245, 55-60.
3. Deisseroth, A., and Dounce, A.L. Catalase: Physical and chemical properties, mechanism of catalysis, and physiological role. *Physiol. Rev.*, 50, 319-375 (1970).
4. Zámocký, M., and Koller, F. (1999) *Prog Biophys Mol Biol.* 72, 19-66.
5. Zhou, M., *et al.* (1997) *Anal Biochem.* 253, 162-168.
6. Mohanty, J. G., *et al.* (1997) *J. Immunol Methods* 202, 133-141.
7. Votyakova, T. V., and Reynolds, I. J. (2004) *Arch Biochem and Biophys.* 431, 138-144.
8. Bai, J., *et al.* (1999) *J Biol Chem.* 274, 26217-26224.
9. Tada-Oikawa, S. *et al.* (1999) *FEBS Lett.*, 442, 65-69.
10. Hampton, M.B., and Orrenius, S. (1997) *FEBS Lett.*, 414, 552-556.
11. Kowaltowski, A.J. *et al.* (2000) *FEBS Lett.*, 473, 177-182.
12. Tome, M.E. *et al.* (2001) *Cancer Res.*, 61, 2766-2733.



MSDS (Material Safety Data Sheet) available online

## Limited Warranty

Assay Designs, Inc. warrants that at the time of shipment this product is free from defects in materials and workmanship. This warranty is in lieu of any other warranty expressed or implied, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose.

Assay Designs must be notified of any breach of this warranty within 48 hours of receipt of the product. No claim shall be honored if Assay Designs is not notified within this time period, or if the product has been stored in any way other than outlined in this publication. The sole and exclusive remedy of the customer for any liability based upon this warranty is limited to the replacement of the product, or refund of the invoice price of the goods.

Assay Designs, Inc. tel: 734.668.6113 or 800.833.8651  
5777 Hines Drive fax: 734.668.2793  
Ann Arbor, Michigan info@assaydesigns.com  
USA www.assaydesigns.com

technical support: technical@assaydesigns.com



## Kits for Charity™

Assay Designs has always been an active contributor to a number of local, national and international charities. We have broadened our charitable contributions by implementing a program called Kits for Charity™.

Each quarter, Assay Designs will feature a different non-religious and non-political charitable organization on our website. For each kit sold during this time period, we will make a monetary contribution to the featured charity. Please check our website for the current quarter's charity to see what organization your purchases are helping to support.

If you have any suggestions for future Kits for Charity™ recipients, please contact us at 800.833.8651 or 734.668.6113.

## Contact Us

For more details concerning the information within this kit insert, or to order any of the Assay Designs' products, please call (734) 668-6113 between 8:30 a.m. and 5:30 p.m. EST. Orders or technical questions can also be transmitted by fax or e-mail 24 hours a day.

Catalog No. 25-0607



Last Revised: 28 June 2007

© 2007

